

WE MADE EXTRAORDINARY
P O W E R



KARINCA 208
Surgical Motor Systems
Battery Operated

AYGUN®
SURGICAL INSTRUMENTS CO., INC.

" Meeting point for surgery..."

EDITON 01



“ Meeting point for surgery...”

Aygün has the right to modify and develop the products declared in this catalogue and its production program.

The weight and dimensional information given in our catalog are approximate values.

KA-4023-Y3 Copyright ©
Aygün Co., Inc 2020 All Right
Reserved, Printed in Türkiye

www.aygun.com

OUR HISTORY

AYGÜN, one of the leading establishments of Turkey in the field of Medical Industry, was founded in 1988. It manufactures Container Systems for sterilization, Surgical Motor Systems, Medical Illumination Systems and Single Use Medical Systems in its facility at Samsun Central Organized Industrial Zone.

AYGÜN meets the needs of significant part of Turkey as a result of its continuous R&D works and the capacity increasing. It manufactures its products in international standards and put them to the world market with a high competitive power. Moreover with its integrated quality management system, it proves that it is a specialized establishment.

Besides being a manufacturer, AYGÜN is also a designer company and with its continuous product development works it increases number of its patents, on the other hand, is contributing career development of rising generation by acting as a school.





STERILIZATION CONTAINER SYSTEMS

All AYGÜN sterilization container systems are designed for perfect steam penetration, traceability, safety and stability. The locking mechanisms secure excellent grip and ensure up and down movement of the lids under steam pressure during the sterilization. We service our bottoms and lids made of 2 mm high quality special alloy aluminum raw material to the whole world with 8 basic color options.

CONTAINER MODELS

STANDARD MODEL
BIO-BARRIER MODEL
PPSU MODEL
WINDOW MODEL
B3 MODEL
PLASMA MODEL
A1 MODEL
SAFE MODEL



SURGICAL MOTOR SYSTEMS

The expertise of AYGÜN in surgical motor systems gathered up many designs and systems that have special structures servicing from neuro and spinal surgery to many fields like cardiovascular, orthopedic, plastic, hand & foot and maxillofacial surgery. The quality, power and ergonomic design of the new generation of high-tech battery powered surgical motor systems once again demonstrates AYGÜN's expertise in the medical field.

SURGICAL MOTOR SYSTEM MODELS

KARINCA 206
KARINCA 208
DERMATOME BES
DERMATOME ST-5
SKIN GRAFT MESHER
PRAXIS 3
PRAXIS 3B
Sternal Saw 3B
HSM Micro System
HSM 2 Neuro & Spine System
Mini Bone 40.000 rpm
Mini Bone Monoblock
Arthroscopy Shaver
Mini Surgical Motor
Mini Compact
Neuro Air Power
Standard Macro
EPS Motor System
Macro Compact
Mobil Stator



SURGICAL ILLUMINATION SYSTEMS

AYGÜN is furthering a number of R&D projects simultaneously thanks to its advanced R&D infrastructure, team of engineers and technical personnel specialized in their field. Thus, the company holds tens of patents and models actively used in the medical sector.

ILLUMINATION SYSTEMS FEATURES

I-LIGHT LD Specifications
Superior Optical performance
Light beams configuration to provide shadowless illumination
Almost perfect color rendition (Typ. 95)
TFT touch screen for lighting control
Adjustable beam sizes to user preferences
Suitable for use with Laminar Flow functional ceilings
Easy to clean
Minimum Energy Consumption
60.000 hours LED life
Indirect Endo light for endoscopic operations



DISPOSABLE STERILIZATION PRODUCTS

AYGÜN continues expanding its product range by adding other single use medical products to the completely silicone medical drainage systems needed by surgery. With its ISO 7 clean room and by carrying its expertise on sterilization and operation solutions one step further, AYGÜN offers professional service to its valued customers in disposable medical items as well.

DISPOSABLE PRODUCTS MODELS

Reusable Laryngeal Mask
Silicone Penrose Drain Tube
Spiral In Tubes
Silicone Nutrition Tube
Silicone Round Drain
Silicone Flat-Soft Drain
Silicone Flat Drain
Intubation Tube
PVC Laryngeal Mask
Silicone Penrose Drain
Silicone Thorax Drain
Tracheostomy tube
Intubation Tube
Vac Silicone Drain
Oxygen Mask

GROWING UP TOGETHER

AYGÜN has an actively growing position in the sector through using a high technology and innovative institutional structure.



OUR VISION

The vision of AYGÜN is to grow in health sector according to the developing technology and sectoral innovations and serve humanity by producing new project and products appropriate to these innovations.

OUR MISSION

AYGÜN builds projects that can contribute to human health all over the world and integrates these projects to its production area and products in the light of technological developments. It constantly stays active with its innovative and explorer soul and remains taking part among the leader companies ruling the sector.

QUALITY MANAGEMENT SYSTEM

- ISO 9001:2015
- ISO 14001:2015
- OHSAS 18001:2007
- ISO 13485:2016
- 93/42/EEC
- FDA 510 (K)





AYGÜN®
SURGICAL INSTRUMENTS CO., INC.


PRODUCTION

Sense of high-grade quality in each process from raw material supply to the shipment is combined with integrated quality management system according to the standards. A big stuff gathered under the same roof and adopted importance of teamwork carry all stages above the dreams and makes what is desired real.

TECHNICAL INFRASTRUCTURE

CNC Machines
Hydraulic and Eccentric Presses
Hot Presses
Conventional CNC Machines
Mold Design and Molding Technology
3D Modelling Technology
Thermoforming
Cutting Machines

Electronic Design and Manufacturing Engineering
Software Technology
Fine Mechanic Techniques
Hardening
Laser Technology
Silicone Extrusion Technology
Clean Room – ISO 7
Chemical Coating and Anodizing Technology












COMBINATION OF PRECIOUS HIGH TECHNOLOGY PRODUCTION MACHINES AND PERFECT CRAFTSMANSHIP










New product development desire and need began in minds becomes sharable and interpretable with 3D printers and technical drawings. After setting the models and testing the prototypes, the production process starts for the approved designs. The production starts with combination of high technology production machines and perfect craftsmanship.

RESEARCH, DEVELOPMENT, QUALITY

AYGÜN increases its power with its endless R&D works. It makes all needed investments so as to research and develop even the smallest part of the complicated products manufactured and carries out successful projects with its professional engineer team. It both develops the existing products with today's technology and performs serious works on new product designs.


TABLE OF CONTENTS

	17 DRILL ATTACHMENTS
	18-19 SHANK DRILL BIT ATTACHMENTS
	20-21 REAMER ATTACHMENT
	22-23 ACETABULAR REAMER ATTACHMENT
	24 QUICK - ACTION CHUCK
	25 RECIPROCATING SAW ATTACHMENT
	26-29 SAGITTAL SAW ATTACHMENT
	30 STERNUM SAW ATTACHMENT
	31 OSCILLATING SAW ATTACHMENT

32 CRANIAL PERFORATOR ATTACHMENT	
33 CRANIAL BURR ATTACHMENT	
34 ORTHOPEDIC SET	
35 DRILL, SAGITTAL SET	
36 DRILL SET	
37 SAGITTAL SAW SET	
38 RECIPROCATING SAW SET	
39 STERNAL SAW SET	
40 CARDIO VASCULAR SET	

Power, Precision and Innovation...


- Less Weight, High Power
- Next Generation Cutting Technology
- New Generation Surgical Motor Technology
- Li-ion Battery Technology



Karınca 208 power box records all data about usage for technical service information. Moreover the system provides users to observe battery level without removing the power box. Karınca 208 can also provide information about battery process cycles like how many times the battery is half charged or fully charged and total operation time in hours, minutes, seconds as well.



Power box oled screen enables users to indicate battery capacity by pressing the test button.



Battery, motor and all electronic components are combined together and exempt from sterilization by aim of increasing life time of the device

KARINCA 208

- Light weight Karınca 208 provides efficient and reliable cutting performance with new generation powerful motor;
- One handpiece solution for reaming, high speed drilling with high torque performance;
- Lithium-ion battery technology provides consistent power during the operation;
- Easy handling and ergonomic design of Karınca 208 ensures reliable and better controlled operation results;
- Efficient torque for acetabular reaming;
- Karınca 208 is an innovative new generation power tool with high quality engineering materials;
- Providing efficiency with engineering solutions;
- Extended battery life;
- Comfortable ergonomic design;
- Easy attachment assembly;
- High speed drilling and oscillating.



"Mechanical structure of Karınca 208 Handpieces provides users non-restrictive washing and disinfection."

BATTERY CHARGING UNITS



- M208-410**
Power box with battery pack
- Powerful brushless DC motor technology
 - Battery capacity: 2.5 Ah / 37 Wh
 - Voltage: 14.8 Volts
 - Integrated oled screen and audible alert notification
 - Not autoclavable
 - 75 x 61 x 160 mm / 540 g



- M206-212**
Autoclavable sterile funnel for aseptic battery loading
- Steam sterilisable
 - 30 x 74 x 104 mm / 55 g



Sterile funnel



- M206-411**
Battery pack, spare part (for service needs only)
- Battery capacity: 2.5 Ah / 37 Wh
 - Voltage: 14.8 Volts
 - Not autoclavable
 - 68 x 56 x 35 mm / 240 g



- M206-226**
Charger unit for 1 battery
- Main voltage 100-240 Volts, 50-60 Hz
 - Charging current 2 A
 - Charging voltage 16.8 Volts
 - Protections: Short circuit / Overload / Over voltage
 - Charging indicator led
 - Battery not included
 - 130 x 250 x 130 mm / 1300 g



- M206-227**
Charger unit for 2 batteries
- Main voltage 100-240 Volts, 50-60 Hz
 - Charging current 2 A (x2)
 - Charging voltage 16.8 Volts
 - Protections: Short circuit / Overload / Over voltage
 - Charging indicator led
 - Battery not included
 - 230 x 250 x 130 mm / 2150 g




- M206-229**
Charger unit for 4 batteries
- Main voltage 100-240 Volts, 50-60 Hz
 - Charging current 2 A (x4)
 - Charging voltage 16.8 Volts
 - Protections: Short circuit / Overload / Over voltage
 - Charging indicator led
 - Battery not included
 - 430 x 250 x 130 mm / 3850 g

Power, Precision and Innovation...

HANDPIECE



M208-200
Hexa - Function Handpiece
Dual trigger system for speed, direction
and oscillation control.
Steam sterilisable
Cannulated Ø 4 mm
77 x 133 x 215 mm / 990 g



Drill / reamer and sagittal / sternum
options in one-handpiece



Forward, reverse and oscillation
mode options provide versatile
alternatives to the users.



Low / high speed options for both
sagittal / sternum modes



DRILL ATTACHMENTS

M208-308
Drill attachment with keyed chuck
Ø4 mm Cannulated
Ø 1.0 mm - 4.0 mm
900 rpm / 4 Nm
30 x 113 mm / 325 g



M8-473
Chuck key
for M208-308

M208-309
Drill attachment with keyed chuck
Ø4 mm Cannulated
Ø 1.0 mm - 6.35 mm
900 rpm / 4 Nm
30 x 126 mm / 370 g




M7-250
Chuck key
for M208-309

M208-313
Drill attachment with hybrid chuck
Both keyed and keyless operation
Ø4 mm Cannulated
Ø 0.3 mm - 3 mm
900 rpm / 4 Nm
30 x 130 mm / 395 g



M7-253
Chuck key
for M208-313


M208-310
Drill attachment with hybrid chuck
Both keyed and keyless operation
Ø4 mm Cannulated
Ø 1.0 mm - 7.4 mm
900 rpm / 4 Nm
34 x 142 mm / 485 g



M7-252
Chuck key
for M208-310

DRILL BITS

Code	(Ø) Diameter mm	Length mm
M7-500	1.5	60
M7-501	2.0	100
M7-502	2.5	100
M7-503	2.7	100
M7-504	3.0	100
M7-505	3.2	100
M7-506	3.5	100
M7-507	4.0	110
M7-508	4.5	130
M7-509	5.0	160
M7-510	6.0	160




Power, Precision and Innovation...



SHANK DRILL BIT ATTACHMENTS

M208-360
Attachment for SQUARE shank
900 rpm / 4 Nm
30 x 114 mm / 270 g



M208-361
Attachment for AO shank
900 rpm / 4 Nm
30 x 106 mm / 260 g



M208-362
Attachment for TRINKE shank
900 rpm / 4 Nm
30 x 104 mm / 270 g





SHANK DRILL BITS

Length mm	(Ø) Diameter mm	SQUARE	AO	TRINKE
38	1.5	M7-520	M7-540	M7-560
45	2.0	M7-521	M7-541	M7-561
53	2.5	M7-522	M7-542	M7-562
57	2.7	M7-523	M7-543	M7-563
57	3.0	M7-524	M7-544	M7-564
57	3.2	M7-525	M7-545	M7-565
65	3.5	M7-526	M7-546	M7-566
70	4.0	M7-527	M7-547	M7-567
70	4.5	M7-528	M7-548	M7-568
95	5.0	M7-529	M7-549	M7-569
120	6.0	M7-530	M7-550	M7-570


M103-375
Attachment for OSTEO shank
Ø4 mm Cannulated
175 rpm / 20 Nm
30 x 89 mm / 225 g





M103-376
Attachment for DHS shank
Ø4 mm Cannulated
175 rpm / 20 Nm
30 x 97 mm / 220 g





M103-377
Attachment for AO-PROTEK shank
Ø4 mm Cannulated
175 rpm / 20 Nm
30 x 85 mm / 260 g



M103-378
Attachment for HUDSON shank
Ø4 mm Cannulated
175 rpm / 20 Nm
30 x 85 mm / 230 g





M103-379
Attachment for HARRIS shank
Ø4 mm Cannulated
175 rpm / 20 Nm
30 x 90 mm / 227 g



M103-380
Attachment for ZIMMER-HUDSON shank
Ø4 mm Cannulated
175 rpm / 20 Nm
30 x 90 mm / 235 g



M103-381
Attachment for ZIMMER shank
Ø4 mm Cannulated
175 rpm / 20 Nm
30 x 90 mm / 235 g





REAMER ATTACHMENT

M103-330
Reamer attachment for wire and pin
Cannulated Ø 4 mm
175 rpm / 20 Nm
30 x 104 mm / 310 g



M7-905
Adaptor for AO-Protek shank
Cannulated Ø 4 mm




M7-908
Adaptor for Hudson
Cannulated Ø 4 mm



M7-907
Adaptor for Harris
Cannulated Ø 4 mm




M7-906
Universal three-jaw chuck
with key M7-250
Max Ø 6.35 mm
Cannulated Ø 4 mm




REAMERS

Intramedullary reamers, flexible, with drilling heads		Guides, blunt			Guides, with drill tip		
Ø mm	Length 430 mm	Ø mm	Length 600 mm	Length 990 mm	Ø mm	Length 600 mm	Length 990 mm
5.5	M7-920						
6.0	M7-921						
6.5	M7-922	1.9	M7-960	M7-961	1.9	M7-962	
7.0	M7-923						
7.5	M7-924						
8.0	M7-925						
8.5	M7-926						
9.0	M7-927						
9.5	M7-928						
10.0	M7-929						
10.5	M7-930						
11.0	M7-931						
11.5	M7-932						
12.0	M7-933	3.4	M7-963	M7-964	3.4		M7-965
12.5	M7-934						
13.0	M7-935						
13.5	M7-936						
14.0	M7-937						
14.5	M7-938						
15.0	M7-939						
15.5	M7-940						
16.0	M7-941						
16.5	M7-942						
17.0	M7-943						
17.5	M7-944						
18.0	M7-945						
18.5	M7-946						
19.0	M7-947						




ACETABULAR REAMER ATTACHMENT

M103-320
Acetabular reamer attachment
100 rpm / 35 Nm
30 x 120 mm / 380 g



M7-705
Adaptor for AO-Protek shaft




M7-706
Adaptor for DIN 58809 shaft



M7-707
Adaptor for Hudson shaft



M7-708
Universal three-jaw chuck
with key M7-250



ACETABULAR REAMER TOOLS



M7-709
Shaft only,
With Harris Connection
for M7-750, M7-757



M7-789
Movable protector case,
Fits with all reamer shafts

Only Acetabular
Reamer



M7-750	Ø 46 mm
M7-751	Ø 48 mm
M7-752	Ø 50 mm
M7-753	Ø 52 mm
M7-754	Ø 56 mm
M7-755	Ø 60 mm
M7-756	Ø 64 mm
M7-757	Ø 68 mm



QUICK - ACTION CHUCK
(for K-wire operations)




RECIPROCATING SAW ATTACHMENT

M103-349
Quick - action chuck
Cannulated Ø 0.6 mm - 1.0 mm
700 rpm / 4 Nm
30 X 87 mm / 350 g

M103-352
Quick - action chuck
Cannulated Ø 0.7 mm - 2.5 mm
700 rpm / 4 Nm
30 X 87 mm / 350 g

M103-353
Quick - action chuck
Cannulated Ø 1.5 mm - 4.0 mm
700 rpm / 4 Nm
30 X 87 mm / 350 g



M103-235
Reciprocating saw attachment
(Without saw blades)
High : 11.500 cpm / \longleftrightarrow 4 mm
Low : 10.000 cpm / \longleftrightarrow 4 mm
33 x 134 mm / 400 g



RECIPROCATING SAW BLADES

M103-236
Usable length: 62 mm
Width: 8 mm



M103-237
Usable length: 75 mm
Width: 13 mm



M103-238
Usable length: 60 mm
Width: 13 mm



M103-239
Usable length: 75 mm
Width: 13 mm



Power, Precision and Innovation...



SAGITTAL SAW ATTACHMENT

M103-300
Sagittal saw attachment
(Without saw blades)

High : 11.500 cpm / $\angle 4^{\circ} 50'$
Low : 10.000 cpm / $\angle 4^{\circ} 50'$
62 x 111 mm / 460 g




M103-301
Sagittal saw long attachment
(Without saw blades)


High : 11.500 cpm / $\angle 4^{\circ} 50'$
Low : 10.000 cpm / $\angle 4^{\circ} 50'$
62 x 175 mm / 700 g




SAGITTAL SAW BLADES



M103-1110
Material Thickness : 0.4 mm
Cutting Thickness : 0.6 mm




M103-1111
Material Thickness : 0.4 mm
Cutting Thickness : 0.6 mm




M103-1112
Material Thickness : 0.4 mm
Cutting Thickness : 0.6 mm


M103-1113
Material Thickness : 0.4 mm
Cutting Thickness : 0.6 mm




M103-1114
Material Thickness : 0.4 mm
Cutting Thickness : 0.6 mm




M103-1115
Material Thickness : 0.6 mm
Cutting Thickness : 0.8 mm




M103-1116
Material Thickness : 0.6 mm
Cutting Thickness : 0.8 mm



M103-1117
Material Thickness : 0.6 mm
Cutting Thickness : 0.8 mm




M103-1118
Material Thickness : 1.19 mm




Power, Precision and Innovation...


M103-1119
Material Thickness: 1 mm




Material Thickness
M103-1136 1.27 mm
M103-1137 1.37 mm




M103-1120
Material Thickness: 1 mm




M103-1138
Material Thickness: 1.27 mm




Material Thickness
M103-1121 1 mm
M103-1122 1.19 mm
M103-1123 1.27 mm
M103-1124 1.37 mm




Material Thickness
M103-1139 1 mm
M103-1140 1.27 mm




Material Thickness
M103-1125 1.19 mm
M103-1126 1.27 mm
M103-1127 1.37 mm




Material Thickness
M103-1141 1 mm
M103-1142 1.19 mm
M103-1143 1.27 mm
M103-1144 1.37 mm
M103-1145 1.47 mm




Material Thickness
M103-1128 1 mm
M103-1129 1.19 mm
M103-1130 1.27 mm
M103-1131 1.37 mm
M103-1132 1.47 mm




M103-1146
Material Thickness: 1.27 mm



Material Thickness
M103-1133 1 mm
M103-1134 1.19 mm
M103-1135 1.27 mm



M103-1147
Material Thickness: 1.19 mm





STERNUM SAW ATTACHMENT

M103-240

Sternum saw attachment
(Without saw blade and guide)

High : 11.500 cpm / \leftrightarrow 4 mm
Low : 10.000 cpm / \leftrightarrow 4 mm
33 x 147 mm / 550 g



OSCILLATING SAW ATTACHMENT

(Re-operation)

M103-300


Sagittal saw attachment
(Without saw blades)

High : 11.500 cpm / \angle 4° 50'
Low : 10.000 cpm / \angle 4° 50'
62 x 111 mm / 460 g



M103-715

Adult saw guide



Adult saw blade	Code	Length	Usable Length	Blade Thickness	Cutting Thickness
	M103-709	58 mm	35 mm	0.6 mm	1 mm
	M103-710	58 mm	35 mm	0.8 mm	1.2 mm


M103-716

Pediatric saw guide



Pediatric saw blade	Code	Length	Usable Length	Blade Thickness	Cutting Thickness
	M103-711	40 mm	18 mm	0.6 mm	0.6 mm
	M103-712	40 mm	18 mm	0.6 mm	0.8 mm

BLADES (for M103-300)



Code	Length	Width	Thickness
M103-390	35 mm	25 mm	0.8 mm
M103-391	50 mm	40 mm	0.8 mm
M103-392	35 mm	75 mm	0.8 mm




CRANIAL PERFORATOR ATTACHMENT
(Dura Stop Perforator Attachment)



CRANIAL BURR ATTACHMENT
(Trauma Operations)


M208-250
Cranial perforator attachment
(without perforator)

900 rpm / 4 Nm
30 x 123 mm / 350 g



◀ For Hudson

M103-255
Cranial burr attachment
Dura guard, rotatable
12.500 rpm
35 x 118 mm / 320 g



TREPANS

- Skull perforator
- With step-down drilling
- Automatic dura stop
- No key needed



M7-120 Ø 6/9 mm
23 x 92 mm / 116 g



M7-121 Ø 9/12 mm
23 x 92 mm / 120 g



M7-122 Ø 12/15 mm
23 x 92 mm / 125 g

CRANIAL BURR TOOLS




M103-203
Tapered cutter



M103-332
Dura guard, rotatable

ORTHOPEDIC SET



M208-1001 Orthopedic Set

Set Content	
M208-200	Hexa - Function Handpiece
M208-410	Power box
M206-212	Autoclavable sterile funnel for aseptic battery loading
M206-226	Charger unit for 1 battery
M208-309	Drill attachment with keyed chuck
M7-502	Drill bit
M7-250	Chuck key
M103-300	Sagittal saw attachment
M103-1146	Sagittal saw blade
M103-320	Acetabular reamer attachment
M103-330	Reamer attachment for wire and pin
M103-352	Quick - action chuck
M103-235	Reciprocating saw attachment
M103-239	Reciprocating saw blades
M14-816	Spray oil
K3-2723-S5	Basket with silicone fixations 546 x 259 x 114 mm
K6-022	Standard Model Container 600 x 285 x 155 mm



M208-410
Power box



M206-226
Charger unit for 1 battery



M14-816
Spray oil, for motor and additional pcs.
After cleaning of every use. Oil inside of all moving pcs of device by spraying method and run the device for 2-3 seconds.




K6-022
Standard Model Container
Box (non perforated) Lid perforated
600 x 285 x 155 mm



K3-2723-S5
Basket with silicone fixations
546 x 259 x 114 mm

DRILL, SAGITTAL SET



M208-1003 Drill, Sagittal Set

Set Content	
M208-200	Hexa - Function Handpiece
M208-410	Power box
M206-212	Autoclavable sterile funnel for aseptic battery loading
M206-226	Charger unit for 1 battery
M208-309	Drill attachment with keyed chuck
M7-502	Drill bit
M7-250	Chuck key
M103-300	Sagittal saw attachment
M103-1146	Sagittal saw blade
M14-816	Spray oil
K3-2731-S6	Basket with silicone fixations 411 x 259 x 114 mm
K6-252	Standard Model Container 475 x 285 x 155 mm



M208-410
Power box




M206-226
Charger unit for 1 battery



M14-816
Spray oil, for motor and additional pcs.
After cleaning of every use. Oil inside of all moving pcs of device by spraying method and run the device for 2-3 seconds.




K6-252
Standard Model Container
Box (non perforated) Lid perforated
475 x 285 x 155 mm



K3-2731-S6
Basket with silicone fixations
411 x 259 x 114 mm

DRILL SET



M208-1005 Drill Set

Set Content	
M208-200	Hexa - Function Handpiece
M208-410	Power box
M206-212	Autoclavable sterile funnel for aseptic battery loading
M206-226	Charger unit for 1 battery
M208-309	Drill attachment with keyed chuck
M7-502	Drill bit
M7-504	Drill bit
M7-250	Chuck key
M14-816	Spray oil
K3-2731-S7	Basket with silicone fixations 411 x 259 x 114 mm
K6-252	Standard Model Container 475 x 285 x 155 mm



M208-410
Power box



M206-226
Charger unit for 1 battery



M14-816
Spray oil, for motor and additional pcs.
After cleaning of every use. Oil inside of all moving pcs of device by spraying method and run the device for 2-3 seconds.




K6-252
Standard Model Container
Box (non perforated) Lid perforated
475 x 285 x 155 mm



K3-2731-S7
Basket with silicone fixations
411 x 259 x 114 mm

SAGITTAL SAW SET



M208-1007 Sagittal Saw Set

Set Content	
M208-200	Hexa - Function Handpiece
M208-410	Power box
M206-212	Autoclavable sterile funnel for aseptic battery loading
M206-226	Charger unit for 1 battery
M103-300	Sagittal saw attachment
M103-1117	Sagittal saw blade
M103-1146	Sagittal saw blade
M14-816	Spray oil
K3-2731-S8	Basket with silicone fixations 411 x 259 x 114 mm
K6-252	Standard Model Container 475 x 285 x 155 mm



M208-410
Power box




M206-226
Charger unit for 1 battery



M14-816
Spray oil, for motor and additional pcs.
After cleaning of every use. Oil inside of all moving pcs of device by spraying method and run the device for 2-3 seconds.




K6-252
Standard Model Container
Box (non perforated) Lid perforated
475 x 285 x 155 mm



K3-2731-S8
Basket with silicone fixations
411 x 259 x 114 mm

RECIPROCATING SAW SET



M208-1013 Reciprocating Saw Set

Set Content	
M208-200	Hexa - Function Handpiece
M208-410	Power box
M206-212	Autoclavable sterile funnel for aseptic battery loading
M206-226	Charger unit for 1 battery
M103-235	Reciprocating saw attachment
M103-237	Reciprocating saw blades
M103-239	Reciprocating saw blades
M14-816	Spray oil
K3-2731-S9	Basket with silicone fixations 411 x 259 x 114 mm
K6-252	Standard Model Container 475 x 285 x 155 mm



M208-410
Power box




M206-226
Charger unit for 1 battery



M14-816
Spray oil, for motor and additional pcs.
After cleaning of every use. Oil inside of all moving pcs of device by spraying method and run the device for 2-3 seconds.




K6-252
Standard Model Container
Box (non perforated) Lid perforated
475 x 285 x 155 mm



K3-2731-S9
Basket with silicone fixations
411 x 259 x 114 mm

STERNAL SAW SET



M208-1009 Sternal Saw Set

Set Content	
M208-200	Hexa - Function Handpiece
M208-410	Power box
M206-212	Autoclavable sterile funnel for aseptic battery loading
M206-226	Charger unit for 1 battery
M103-240	Sternum saw attachment
M103-715	Sternum saw guide
M103-710	Sternum saw blade
M14-816	Spray oil
K3-2731-S10	Basket with silicone fixations 411 x 259 x 114 mm
K6-252	Standard Model Container 475 x 285 x 155 mm



M208-410
Power box




M206-226
Charger unit for 1 battery



M14-816
Spray oil, for motor and additional pcs.
After cleaning of every use. Oil inside of all moving pcs of device by spraying method and run the device for 2-3 seconds.




K6-252
Standard Model Container
Box (non perforated) Lid perforated
475 x 285 x 155 mm



K3-2731-S10
Basket with silicone fixations
411 x 259 x 114 mm

CARDIO VASCULAR SET

Note:



M208-1011 Cardio Vascular Set

Set Content	
M208-200	Hexa - Function Handpiece
M208-410	Power box
M206-212	Autoclavable sterile funnel for aseptic battery loading
M206-226	Charger unit for 1 battery
M103-240	Sternum saw attachment
M103-715	Sternum saw guide
M103-710	Sternum saw blade
M103-300	Sagittal saw attachment
M103-392	Re-Operation saw blade
M14-816	Spray oil
K3-2731-S11	Basket with silicone fixations 411 x 259 x 114 mm
K6-252	Standard Model Container 475 x 285 x 155 mm



M208-410
Power box




M206-226
Charger unit for 1 battery



M14-816
Spray oil, for motor and additional pcs.
After cleaning of every use. Oil inside of all moving pcs of device by spraying method and run the device for 2-3 seconds.



K6-252
Standard Model Container
Box (non perforated) Lid perforated
475 x 285 x 155 mm



K3-2731-S11
Basket with silicone fixations
411 x 259 x 114 mm

CONTAINER SET FOR STERILIZATION



K8-123V
Lid, B3 Barrier Model
602 x 285 mm


K6-107
Container, non perforated
600 x 272 x 138 mm



K3-742
Basket made from
stainless steel sheet (with
stand)
540 x 253 x 80 mm



K2-462
Wire basket
540 x 250 x 70 mm




K2-526
Silicone mat
520 x 230 mm



K8-343V
Lid, B3 Barrier Model
477 x 285 mm


K6-307
Container, non perforated
475 x 272 x 138 mm



K3-750
Basket made from
stainless steel sheet (with
stand)
405 x 253 x 80 mm



K2-1502
Wire basket
405 x 250 x 70 mm



K2-524
Silicone mat
380 x 230 mm



K2-674
Class 6 Emulating Indicator Card
1 Pack (250 pcs)




	code	color
Basic Set	K2-700	Red
Thoracotomy	K2-701	Blue
Orthopedic Set	K2-702	Green
Bone Set	K2-703	Yellow
Laminectomy	K2-704	Grey
Neuro Surgery	K2-705	Black
Cardiovascular	K2-706	Pink
Thoracotomy	K2-707	Brown

Identification labels (aluminium)
Coloured identification labels
Engraved to your specifications
Maximum 16 letters



K2-673
Labels (paper)
100 pcs. With indicator.
For small size and special
containers



K2-672
Labels (paper)
100 pcs. With indicator.
For full and half size containers



K1-900
Modular Auto Lock Mechanism (2 Pcs)




K2-660 red
Security seal
100 pcs made of plastic



K2-740
Single use paper
filter 100 pcs
with indicator.
Aygun mark.



For more information and alternative please see
our Sterilization Container System catalog.
From our website www.aygun.com





Standard Model



Bio-Barrier Model



High Performance
Polymer PPSU Lids



B3-PTFE - Valve - Bio-stop



Plasma Model



Window Model



A1 Model



Dental, Mini Model



Mini Bone



Mini Bone
Monoblock



Arthroscopy
Shaver System



Karinca 206




Karinca 208



HSM-2
Neuro System



HSM
Neuro System



HSM
Micro System



Dermatome ST5



Dermatome




Praxis 3B



Praxis 3



Compact 2



Standard Macro



Standard Micro



Electro Pneumatic
Surgical Motor




Macro Compact 2



Mini Compact



Mini Surgical Motor



NAP Neuro Air
Power



INDEX

CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE
K1-900	43	M7-507	17	M7-754	23	M103-235	25	M103-1125	28
K2-462	42	M7-508	17	M7-755	23	M103-236	25	M103-1126	28
K2-524	42	M7-509	17	M7-756	23	M103-237	25	M103-1127	28
K2-526	42	M7-510	17	M7-757	23	M103-238	25	M103-1128	28
K2-660	43	M7-520	18	M7-789	23	M103-239	25	M103-1129	28
K2-672	43	M7-521	18	M7-905	20	M103-240	30	M103-1130	28
K2-673	43	M7-522	18	M7-906	20	M103-255	33	M103-1131	28
K2-674	43	M7-523	18	M7-907	20	M103-300	26, 31	M103-1132	28
K2-700	43	M7-524	18	M7-908	20	M103-301	26	M103-1133	28
K2-701	43	M7-525	18	M7-920	21	M103-320	22	M103-1134	28
K2-702	43	M7-526	18	M7-921	21	M103-330	20	M103-1135	28
K2-703	43	M7-527	18	M7-922	21	M103-332	33	M103-1136	29
K2-704	43	M7-528	18	M7-923	21	M103-349	24	M103-1137	29
K2-705	43	M7-529	18	M7-924	21	M103-352	24	M103-1138	29
K2-706	43	M7-530	18	M7-925	21	M103-353	24	M103-1139	29
K2-707	43	M7-540	18	M7-926	21	M103-375	19	M103-1140	29
K2-740	43	M7-541	18	M7-927	21	M103-376	19	M103-1141	29
K2-1502	42	M7-542	18	M7-928	21	M103-377	19	M103-1142	29
K3-742	42	M7-543	18	M7-929	21	M103-378	19	M103-1143	29
K3-750	42	M7-544	18	M7-930	21	M103-379	19	M103-1144	29
K3-2723-S5	34	M7-545	18	M7-931	21	M103-380	19	M103-1145	29
K3-2731-S6	35	M7-546	18	M7-932	21	M103-381	19	M103-1146	29
K3-2731-S7	36	M7-547	18	M7-933	21	M103-390	31	M103-1147	29
K3-2731-S8	37	M7-548	18	M7-934	21	M103-391	31	M206-212	14
K3-2731-S9	38	M7-549	18	M7-935	21	M103-392	31	M206-226	15
K3-2731-S10	39	M7-550	18	M7-936	21	M103-709	30	M206-227	15
K3-2731-S11	40	M7-560	18	M7-937	21	M103-710	30	M206-229	15
K6-022	34	M7-561	18	M7-938	21	M103-711	30	M206-411	14
K6-107	42	M7-562	18	M7-939	21	M103-712	30	M208-200	16
K6-252	35	M7-563	18	M7-940	21	M103-715	30	M208-250	32
K6-307	42	M7-564	18	M7-941	21	M103-716	30	M208-308	17
K8-123V	42	M7-565	18	M7-942	21	M103-1110	26	M208-309	17
K8-343V	42	M7-566	18	M7-943	21	M103-1111	26	M208-310	17
M7-120	32	M7-567	18	M7-944	21	M103-1112	26	M208-313	17
M7-121	32	M7-568	18	M7-945	21	M103-1113	27	M208-360	18
M7-122	32	M7-569	18	M7-946	21	M103-1114	27	M208-361	18
M7-250	17, 20	M7-570	18	M7-947	21	M103-1115	27	M208-362	18
M7-252	17	M7-705	22	M7-960	21	M103-1116	27	M208-410	14
M7-253	17	M7-706	22	M7-961	21	M103-1117	27	M208-1001	34
M7-500	17	M7-707	22	M7-962	21	M103-1118	27	M208-1003	35
M7-501	17	M7-708	22	M7-963	21	M103-1119	28	M208-1005	36
M7-502	17	M7-709	23	M7-964	21	M103-1120	28	M208-1007	37
M7-503	17	M7-750	23	M7-965	21	M103-1121	28	M208-1009	39
M7-504	17	M7-751	23	M8-473	17	M103-1122	28	M208-1011	40
M7-505	17	M7-752	23	M14-816	34, 35	M103-1123	28	M208-1013	38
M7-506	17	M7-753	23	M103-203	33	M103-1124	28		



FACTORY

AYGÜN CERRAHİ ALETLER A.Ş.
Organize Sanayi Bölgesi
Yaşar Doğu Caddesi, No: 76
Kerimbey OSB, 55300, Tekkeköy
SAMSUN, TÜRKİYE

T : +90 (362) 266 66 80
F : +90 (362) 266 93 30

aygun@aygun.com
www.aygun.com