

Our Priority Is Your Health

Get Your Protection & Stay Healthy

SURGICAL PROTECTIVE PRODUCTS





Drenay has the right to modify and develop the products declared in this catalogue and its production program.

KA-01-Y1 Copyright® Drenay Co., Inc 2020 All Right Reserved, Printed in Türkiye

www.drenay.com

COMPANY

Our History


DRENAY, which was established with more than 30 years of experience in the medical industry and manufactures sterile packaged single use systems for the needs of the medical sector, is an AYGÜN corporation. DRENAY has production facilities at international standards, is a designer company as well as being a manufacturer with its professional management approach and specialist engineers in their fields.

Our Mission

Our company, which has projects that can contribute to human health worldwide, integrates its projects into its production area and products in the light of technological developments. Drenay is being constantly active with its innovative and open for improvement spirit, is among the leading companies that steer the sector.

Our Vision

Grow in the healthcare sector according to the developing technology and sectoral innovations and to serve humanity by producing new projects and products suitable for these innovations.







	PAGE
INTRODUCTION	6
SURGICAL FACE MASK	9
PARTICULE FILTERING MASK	13 16
FACE SHIELD	17 18
SINGLE USE GOWN	19 21

Production

DRENAY continues its activities in the medical sector with ISO 7 (Class 10000) and ISO 8 (Clas 100000) production areas in compliance with international standards and qualified personnel in their fields.

With our superior quality understanding at all stages from raw material entry to shipment, Our main goal is customer satisfaction.









QUALITY POLICY

As DRENAY MEDICAL SYSTEMS CO., INC. in line with our vision, mission, principles and values, we aim to act in line with the legal requirements through our Integrated Management System, not to lose the trust of our customers and contribute to the improvement of our nation with the sustainable development.

About quality, environment, occupational health and safety issues; An Integrated Management System is established and implemented in order to prevent contamination, injury and health deterioration in accordance with a management understanding, risks and environmental impacts of our operations.

In the light of quality policy;

- To provide our products with our fast and reliable service, to protect the environment, to use our resources such as our experienced workforce, technological infrastructure and materials efficiently and effectively and to be one of the leading companies in our sector by continuously improving our system.
- ٠ To establish management systems to ensure improvement in line with DRENAY policies, principles and objectives.
- To increase customer satisfaction by implementing a management system that starts with the research of ٠ customer and market requests and continue with sustaining all processes ranging from design, production and after sales support.

In the light of Occupational Health and Safety policy;

- Protecting and using natural resources in the most efficient way.
- . in which we operate.
- by researching new technologies.
- To contribute to continuous development by creating environmental goals and objectives.

In the light of Occupational Health and Safety policy;

- To provide a healthy and safe working environment, to prevent occupational accidents and occupational diseases.
- By reviewing the working conditions frequently, applying the measures in the disruptive issues, providing the necessary measures and creating a working principle. Reviewing the policy according to the changing conditions.
- We are committed to working under these policy titles and provide resources on this path.

Sincerely...

٠



To comply with the national and international legal regulations and obligations related to the environment

To take into account the environmental factors in the evaluation of manufacturing and operation processes



HOW TO USE SURGICAL MASK



WASH YOUR HANDS BEFORE WEARING A MASK



LOCATE THE METALLIC STRIP AND PLACE IT ON THE NOSE BRIDGE


PRESS THE METALLIC STRIP

TO FIT THE SHAPE OF THE NOSE

REMOVE THE MASK FROM BEHIND BY

HOLDING THE STRINGS WITH


CLEAN HANDS



CHECK THE NEW MASK TO MAKE SURE IT'S NOT DAMAGED



SECURE THE STRINGS BEHIND YOUR HEAD OR OVER YOUR EARS





REPLACE THE MASK IF IT GETS DAMP AND DO NOTREUSE IT



DISPOSE THE MASK IN A CLOSED BIN WITHOUT TOUCHING THE FRONT



ENSURE THE PROPER SIDE OF THE MASK FACES OUTWARDS




COVER MOUTH AND NOSE FULLY, MAKING SURE THERE ARE NO GAPS



DO NOT TOUCH THE MASK WHILE USING IT, IF YOU DO WASH YOUR HANDS



LEARN HOW TO USE THE MASK PROPERLY IN USER MANUAL






SURGICAL FACE MASK Surgical Protective Products


SURGICAL FACE MASK **SPECIFICATIONS**

- Single use only.
- 3 Ply.
- Non pilling.
- Manufactured from non-woven materials with ultrasonic welding technology.
- EN 14683 Type II, high bacterial filtration efficiency BFE 98%.
- Ergonomic structure fully compatible with the face.
- Comfortable to use thanks to its ideal air permeability.
- White, green and blue color options.
- With earloop or ties.
- Hygienic and hypoallergenic.
- Do not contain latex and fiberglass.

EN 14683 + AC 2019			
Microbial Cleanliness Level (cfu/g)	≤ 30	≤ 30	≤ 30
Splash Resistance (kPa)	Not Required	Not Required	≥ 16,0
Pressure Difference (Pa/cm²)	< 40	< 40	< 60
Bacterial Filtration Efficiency (BFE), (%)	≥ 95	≥ 98	≥ 98
TEST	ΤΥΡΕ Ι	TYPE II	TYPE IIR







1-3	2	4
Spunbond Non Woven	MeltBlown/Spunbond Non Woven	Nose Wire
Material that are used inner and outer mask layer provides water repellency and extra protection.	Material that are used in middle layer of the mask provides bacterial filtration.	Polypropylene coated aluminum wire





REF NO	TYPE	PIECE	COLOR
CM1-L	Type I	1	
CM2-L	Type I	1	
CM3-L	Type I	1	
CMM1-L	Type II	1	
CMM2-L	Type II	1	
CMM3-L	Type II	1	
CMM1-LR	Type IIR	1	
CMM2-LR	Type IIR	1	
CMM3-LR	Type IIR	1	

DIFOF

Green White Blue



Mask with ties

White Blue Green

REF NO	TYPE	PIECE	COLOR
CM1-B	Type I	1	
CM2-B	Type I	1	
СМЗ-В	Type I	1	
CMM1-B	Type II	1	
CMM2-B	Type II	1	
СММЗ-В	Type II	1	
CMM1-BR	Type IIR	1	
CMM2-BR	Type IIR	1	
CMM3-BR	Type IIR	1	




FILTERING MASK



PARTICULE FILTERING MASK Surgical Protective Products

PARTICLE FILTERING MASK SPECIFICATIONS

- Depend on filtration levels FFP1, FFP2 and FFP3 versions available with and without valve.
- It is made of non-woven materials with low breathing resistance and high filtration efficiency.
- Manufactured from non-woven materials with ultrasonic welding technology.
- Ergonomic design provides a perfect fit to the face.
- Do not contain latex and fiberglass.
- Hygienic and hypoallergenic.



1-3-5 Spunbond Non Woven


Provides extra protection and water repellency.

2-4 MeltBlown Non Woven

Material that are used in middle layers of the mask provides bacterial filtration.



Particle Filtration Efficency for 0,3 Microns






REF NO	PIECE	VALVE
FM1-25FFP1	25	-
FM1-50FFP1	50	-
FM1-25FFP2	25	-
FM1-50FFP2	50	-
FM1-25FFP3	25	-
FM1-50FFP3	50	-
FMV1-25FFP1	25	V
FMV1-50FFP1	50	V
FMV1-25FFP2	25	V
FMV1-50FFP2	50	V
FMV1-25FFP3	25	V
FMV1-50FFP3	50	V





DRENAY, which was established with more than 30 years of experience in the medical industry and manufactures sterile packaged single use systems for the needs of the medical sector, is an AYGÜN corporation.





FACE SHIELD SURGICAL PROTECTIVE PRODUCT



FACE SHIELD **SPECIFICATIONS**

- Adjustable, ergonomic design.
- Provides long-term use thanks to the highly transparent and antifogging visor material.
- Protects face mask and face to body fluids like blood and other environmental effects.
- Suitable for use with masks and glasses.
- Provides ease of use thanks to the visor that can be opened and closed. •
- Reusable thanks to replacable visor.





Antifogging Protects from Can be used Extra clear Protects from Protects from with mask pathogens liquids visor particles



SINGLE USE GOWNS SPECIFICATIONS

- It is made of laminated non-woven materials by ultrasonic welding.
- Antistatic and resistant to liquid splashes.
- 4 different size options and inner lining for easy dressing.
- Adhesive zipper closure and elasticized sleeve, leg cuffs, hood and waist area provide extra protection.


REF NO	PIECES	SIZE	STERILE
T1-2536	25	S	
T1-2538	25	М	
T1-2540	25	L	
T1-2542	25	XL	
T1-2544	25	XXL	
T1-25S36	25	S	S
T1-25S38	25	М	S
T1-25S40	25	L	S
T1-25S42	25	XL	S
T1-25S44	25	XXL	S







KEEP YOUR SAFEDISTANCE





FACTORY

Organize Sanayi Bölgesi Yaşar Doğu Caddesi No:76/2 Kerimbey OSB 55300 Tekkekoy, Samsun, Turkey

T : +90 (362) 266 88 44 F : +90 (362) 266 88 43

www.drenay.com