

Power System



MACRO COMPACT Pneumatic Surgical Motor

Edition 6 -

www.aygun.com



Organize Sanayi Bölgesi Yaşar Doğu Cad. No:76 55330 Kutlukent-SAMSUN TÜRKİYE

Tel	: 90 (362) 266 66 80
Fax	: 90 (362) 266 93 30
E-mail	: aygun@aygun.com
Web	:www.aygun.com

 According to the production programme Aygün's rights is reserved in production to improve, modify or delete the products declared in this catalog.

KA - 4011-Y1 Copyright © Aygün Co., Inc 2013 All Right Reserved, Printed in Türkiye



AYGÜN, one of the leading establishment of Turkey in the field of Medical Industry, was founded in 1988.

It manufactures Container Systems for Sterilization, Surgical Motor Systems, Medical Illumination Systems and Single Use Medical Systems in its facility at Samsun Central Organized Industrial Zone.

AYGÜN meets the needs of significant part of Turkey as a result of its continuous R&D works and the capacity increasing. It manufactures its products in international standards and put them to the world market with a high competitive power. Moreover with its integrated quality management system, it proves that it is a specialized establishment.

Besides being a manufacturer, AYGÜN is also a designer company and with its continuous product development works it increases number of its patents, on the other hand, is contributing career development of rising generation by acting as a school.



<complex-block><complex-block><complex-block>


AYGUN

AYGUN




www.aygun.com











M8-510 Spray oil, for motor and components before sterilization.

After cleaning of every use. Oil inside of all moving parts of device by spraying method and run the device for 2-3 seconds. Practice spraying inside of the device until the black used oil comes out and new oil clearly appears out of device as clear color and then put the device in store or process sterilization.





M7-220 3 m M7-221 5 m Combine hoses for air supply (for M7-230 and M7-150)

Sterilization= 3 Atm 134° C



Air pressure reducer, two way

For security; the systems must be operated by using compressed dry air or nitrogen gas, the oxygen gas must never be used.



EPS Motor System - 2 is used when the regulator at the output of pressure welding is directly connected to the handle of function tool with compressed air hose. EPS-2 is a pneumatic surgical motor system with trigger mechanism. The system does not require an additional foot pedal as the speed control is performed with the speed control switch on the handle of the product In the design of EPS Motor System - 2, also a safety lock system has been

developed on speed control trigger to prevent it from being pushed involuntarily.

M7-220

M7-221

PNEUMATIC - SURGICAL

M7-230

MACRO COMPACT Pneumatic Surgical Motor



M4-219



M4-220



M4-221 Adaptor for SQUARE s



M4-222 Adaptor for AO shank




M4-223 Adaptor for TRINKE sha



M4-224 Adaptor for OSTEO sha





2



Attachment with small three - jaw chuck, with key, Ø 4 mm, exchangeable, use with M4-200 and M4-1100


Attachment with three - jaw chuck, with key, Ø 6.35 mm, exchangeable, use with M4-200 and M4-1100

hank		TORS
		CHUCKS AND ADAPTORS
ank	•	CHUCKS /
ank		
14-225 Quick-action chuck fo se with M4-200 and	r boring wires up to Ø 2.5 mm, M4-1100	






EPS Motor System - 2 is used when the regulator at the output of pressure welding is directly connected to the handle of function tool with compressed air hose. EPS-2 is a pneumatic surgical motor system with trigger mechanism. The system does not require an additional foot pedal as the speed control is performed with the speed control switch on the handle of the product In the design of EPS Motor System - 2, also a safety lock system has been developed on speed control trigger to prevent it from being pushed involuntarily.



MACRO COMPACT Pneumatic Surgical Motor



For M4-1200

Code	Working Length	Width
M7-1110	25	5
M7-1111	25	10
M7-1112	35	5
M7-1113	35	10
M7-1114	50	20
M7-1115	50	10
M7-1116	50	15
M7-1117	50	20
M7-1118	50	25
M7-1119	50	30
M7-1120	65	20
M7-1121	65	25
M7-1122	65	30
M7-1123	65	40
M7-1124	65	50
M7-1125	35	75
M7-1126	50	40
M7-1127	75	20
M7-1128	75	30
M7-1129	90	15
M7-1130	90	20
M7-1131	90	30
M7-1132	90	25
M7-1133	90	25
M7-1134	90	25
M7-1135	90	25
M7-1136	105	19




Total Length	Thickness	
50	0.5	
50	0.5	
60	0.5	
60	0.5	
60	0.5	
75	0.5	
75	0.5	
75	0.7	
75	0.7	
75	0.7	
90	0.7	
90	0.7	
90	0.7	
90	0.7	
90	0.7	
60	0.5	
75	0.5	
105	1	
105	1	
114	1	
114	1	
114	1	
114	0.9	
114	1	
114	1.20	
114	1.5	
129	1.20	

SAGITTAL SAW BLADES

5





MACRO COMPACT Pneumatic Surgical Motor

	Intramed flexible, v	ullary reamers, vith drilling heads		Guides, blu	nt		Guides, w	ith drill tip	
	Ømm	lenght 430 mm	Ømm	lenght 600 mm	lenght 990 mm	Ømm	lenght 60 mm	lenght 990 mm	
	5.5 6.0 6.5 7.0 7.5	M7-920 M7-921 M7-922 M7-923 M7-924	1.9	M7-960	M7-961	1.9	M7-962		
	8.0 8.5 9.0 9.5 10.0 10.5 11.0 11.5 12.0 12.5 13.0 13.5 14.0 14.5 15.0 15.5 16.0 16.5 17.0 17.5 18.0 18.5	M7-925 M7-926 M7-927 M7-928 M7-929 M7-930 M7-930 M7-931 M7-932 M7-933 M7-933 M7-934 M7-935 M7-936 M7-937 M7-938 M7-937 M7-938 M7-939 M7-940 M7-941 M7-942 M7-943 M7-944 M7-945 M7-946	3.4	M7-963	M7-964	3.4		M7-965	
 DIN 58809	19.0	M7-947							



Only Acetabular Reamer





/17-750	Ø 46 mm
/17-751	Ø 48 mm
/17-752	Ø 50 mm
/17-753	Ø 52 mm
/17-754	Ø 56 mm
/17-755	Ø 60 mm
/17-756	Ø 64 mm
/17-757	Ø 68 mm

M7-710

Shaft only, With Harris connection For M7-750, M7-757

M7-791 Movable protector case, fitting all reamer shafts sterilization: 3 atm 143° c



7



EPS Motor System - 2 is used when the regulator at the output of pressure

welding is directly connected to the handle of function tool with compressed

air hose. EPS-2 is a pneumatic surgical motor system with trigger mechanism.

The system does not require an additional foot pedal as the speed control

is performed with the speed control switch on the handle of the product In

the design of EPS Motor System - 2, also a safety lock system has been

developed on speed control trigger to prevent it from being pushed involuntarily.

PNEUMATIC - SURGICAL

M7-220


M7-221





M7-230

EPS Motor System - 2 is used when the regulator at the output of pressure welding is directly connected to the handle of function tool with compressed air hose. EPS-2 is a pneumatic surgical motor system with trigger mechanism. The system does not require an additional foot pedal as the speed control is performed with the speed control switch on the handle of the product In the design of EPS Motor System - 2, also a safety lock system has been developed on speed control trigger to prevent it from being pushed involuntarily.



M7-230

EPS MOTOR SYSTEM 2 Sternum reciprocating saw With saw guide (M4-715) and blade (M4-710) No key needed Saw blade M4-710 and Saw guide M4-715 can be mounted in 2 differend positions. 220 x 165 mm / 1440 g

M4-1500

8



PNEUMATIC - SURGICAL







Container Set For Sterilization








K6-307 Container, non perforated 475 x 272 x 138 mm













code
K2-700
K2-701
K2-702
K2-703
K2-704
K2.705

Identification labels (aluminium) Coloured identification labels Engraved to your specifications Maximum 16 letters



K2-673 Labels (paper) 100 pcs. With indicator. For small size and special containers

6	1	STR NO	
	140-1	STR DATE	
BUM	0 10 I	DUE DATE	
1	8°8	NAME	
U	-	STEAM BLACK	2



K1-900



K2-660 red Security seal 100 pcs made of plastic

For more information and alternative please see our Sterile Container System catalog. From our website www.aygun.com






Edition 71

11





Edition 6

Meeting Point For Surgery™



Organize Sanayi Bölgesi Yaşar Doğu Cad. No:76 55330 Kutlukent - SAMSUN / TÜRKİYE

 90 (362) 266 66 80

 E-mail: aygun@aygun.com

1 90 (362) 266 93 30 **Web**: www.aygun.com